emmi[®]PET

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Redaktioneller Stand: November 2022 LK

emmi[®]PET

PROFESSIONAL DENTAL CLEANING FOR ANIMALS WITH 100% ULTRASOUND

User Manual - emmi[®]-pet

Thank you for choosing our emmi[®]-pet ultrasonic toothbrush for dental cleaning and oral hygiene for animals. To get familiar how to use our ultrasonic device for dental cleaning, emmi[®]-pet, please read this user manual carefully and completely. This secures you reach the best possible results for your animal.

Content emmi[®]-pet:

- 1x emmi[®]-pet Handpiece,
- 1x emmi[®]-pet 100-240-V Charging station,
- 1x emmi[®]-pet A1/A2 Ultrasound attachments,
- 1x emmi[®]-pet ZA Ultrasound toothpaste for animals (75ml).

CORRECT POSITIONING

APPLICATION VIDEO

Premium Partner des VDH

HEALTHY TEETH - HEALTHY ANIMAL

Important - Charging of batteries / Performance per charge

Please make sure that the handpiece is switched off during charging, otherwise charging will not be possible. If you use the unit 10 minutes once a day, you can use it on the battery for approx. 10 days. Afterwards, charging as needed is required.

First Use

1. Always insert the ultrasound attachment in vertical direction and remove it also vertically. This avoids damage to the handpiece contact pin due to jamming.
2. Connect the plug to the socket, place the handpiece with ultrasound attachment on the charging station and charge the unit for 24 hours.

LED 1 lights up blue as soon as the brush is placed on the charging station. Bright blue light means: charging in progress, dim blue light means: charging complete.

IMPORTANT SECURITY ADVISE

The ultrasonic devices of Emmi Ultrasonic GmbH reflect the newest and most effective way of dental cleaning and oral hygiene for humans and animals. Please read the user manual carefully and completely as well as the enclosed warning notices. To avoid an electric shock please pay attention to the below:

1. Always place the charging station in such a way that it cannot fall into bathtubs or washbasins
2. Use the charging station on dry subsoil and in dry surrounding only.
3. In case the charging station has fallen into water, don't touch it.
4. Interrupt immediately the power supply.
5. Never use the charging station in the bath tub.

LED Display

LED 1: Red, green and blue light on the button in one LED.

1. LED blue – CHARGING

LED 1 lights up blue as soon as the brush is placed on the charging station. Bright blue light means: charging in progress, dim blue light means: charging complete.

2. LED green – ON

If the unit is switched on and the ultrasound attachment is inserted, the green LED is lit continuously.

3. LED red blinking - CHECK / FAULT / BATTERY CHARGE

If there is a fault, the red LED lights up and the device switches off. If the battery is heavily discharged, LED 1 flashes red a few times after the cleaning process. Now the brush should be charged without fail, as the battery energy may no longer be sufficient for another cleaning process.

4. Unit switched off - No display

Additional features

- Attachment detection (Transport Safe function): without an attachment, the unit switches off automatically and cannot be switched on. The green LED lights up briefly and switches off again.
- Automatic switch-off (Akku Safe function): the device switches off automatically after 6 minutes.

During charging, always insert the ultrasound attachment for a better control of the unit condition.

Intended Use

To avoid the risk of damages, fire, electrocution or injuries please follow the below:

Use this product only for its intended use as described in this manual. This means, use the device only for the dental cleaning and oral hygiene of animals and only in conjunction with other devices and components recommended by Emmi Ultrasonic GmbH. Any other or beyond usage is considered as not intended. In case of not intended use Emmi Ultrasonic GmbH refuses the liability for possibly appearing damages and the guarantee for a functional operation of the components. Intended use also includes following all indications stated in this manual. This is pointed out by the symbol ⚠ on the type plate. A correct and secure operation of the product requires appropriate transport, storage, placement and set up as well as an elaborate handling and maintenance. The recurring testings according to DGUVV3 are to be considered. The device needs to be positioned in a way that the power plug is freely accessible at any time.

Not Intended Use

In case of not intended use dangers for life and body as well as material damages can occur.

For this reason not intended use is permitted. In case of not intended use Emmi Ultrasonic GmbH refuses the liability for possibly appearing damages and the guarantee for a functional operation of the components. To avoid not intended use, please follow the below.

- This device for dental cleaning and oral hygiene may be operated only with the supply voltage printed on the type plate (bottom side of the charging station).
- Don't use any flammable or aggressive chemicals, acids, bases, etc. to clean the device.
- In case transport damages are experienced, don't connect the charging station to the power supply.
- The device may be operated only indoor at room temperature.

HEALTHY TEETH - HEALTHY ANIMAL

One way to get your animal used to dental and oral care with the original Ultrasound of emmi®-pet (you know your pet best after all)

- 1 Press 1 cm of Ultrasound toothpaste on your finger. Let your animal sniff and lick the cream.
- 2 Switch the unit on and show it to your animal.
- 3 Hold the handpiece with the Ultrasound tip near the mouth of your animal and talk reassuringly to your animal.
- 4 Lift carefully his lips and introduce slowly the ultrasonic attachment into his mouth, in order to get your animal used to the hygienic procedure.
- 5 Now start cleaning the teeth. Go from tooth to tooth and repeat the following application steps. With the bristles, distribute loosely the Ultrasound toothpaste over the inner and outer surfaces of teeth and chewing surfaces. Start with his back teeth while holding the bristles of the ultrasonic attachment and tip motionless on the teeth without any pressure. The ultrasonic attachment allows you to cover 1 to 3 teeth at once. Leave the brush approx. 10 seconds on the teeth and move on to the next teeth. Please do not move the unit during cleansing, otherwise it will not function. The micro-bubbles that are generated in the Ultrasound toothpaste by the emmi®-pet Ultrasound ensure a thorough and friction-free cleansing.
- 6 After cleansing, wash the bristles under running water. The unit remains switched on during this procedure. After rinsing, you can switch off the unit. The Ultrasound has a sterilization effect which ensures disinfection after rinsing. After each cleansing, the device can be repositioned on the charging station. This does not damage the battery.

100% ULTRASOUND TO PREVENT DENTAL DISEASES

Daily Use

We recommend to use the emmi[®]-pet for dental and mouth care approximately 1 to 3 times a week (veterinarians recommend a daily use). The duration of the cleaning should be between 3 - 10 minutes depending on the degree of contamination of your pets teeth. If a regular cleaning is happening (daily tooth cleaning) 3 minutes are enough as part of the prevention. If your pets teeth are cleaned irregularly with the emmi[®]-pet, we recommend intensive sonication of the teeth for up to 10 minutes. Ultrasound is absolutely harmless to your animal.

How does emmi[®]-pet work in combination with the special emmi[®]-pet ultrasonic toothpaste enriched with microbubbles?

The emmi[®]-pet Ultrasound toothpaste for animals has been developed specially for dental care in animals. The emmi[®]-pet, along with the Ultrasound unit for dental and mouth care, gently cleans the teeth of your animal without any friction. It is important that you use the emmi[®]-pet Ultrasound toothpaste for animal use for the dental care of your animal. The supplied emmi[®]-pet ultrasonic toothpaste contains microscopic air bubbles (microbubbles) that have been enriched in an additional production step. The Goethe University of Frankfurt am Main has confirmed this. More than 3.5 million microbubbles were detected in a range of one mm³. The ultrasonic vibrations destroy the structure of the microbubbles and they collapse. This creates a suction effect that removes the biofilm and the plaque from the teeth. Afterwards we recommend brushing the bristles over the teeth to remove the loosened biofilm and the plaque from the tooth. All of this happens without your animal noticing it. A traditional toothpaste will not generate the micro-activity (generation of micro-bubbles) which is indispensable for the functioning of the emmi[®]-pet Ultrasound in dental care applications.

IMPORTANT:

The device has a vibration and noiseless functioning in order not to frighten or upset your animal. The correct functioning of the unit is ensured if the green display under the switch of the handpiece is continuously lit during operation. The unit cannot be used if the Ultrasound tip is not inserted.

Acute inflammations of the gums

In case of acute problems, hold the bristles loosely on the concerned teeth and then directly on the gingiva. Leave them there for about 30 to 60 seconds. This allows you a quick and effective improvement. Prolonged use is not critical, since the Ultrasound is absolutely harmless.

Preventing a bad smell from the mouth

Since the emmi®-pet Ultrasound reduces the bacteria on and in the gingiva, it also prevents the bad smell caused by living bacteria. It is effective for multiple hours.

How the Emmi®-pet Ultrasound works

The dental and mouth care with Ultrasound is an efficient household method. The emmi®-pet functions differently from the traditional, manual, electric or sound toothbrushes using mechanical friction. The emmi®-pet applies the capabilities of modern Ultrasound technology. During the dental care with emmi®-pet, the cleansing is not produced by the usual «grinding» of traditional toothbrushes, but by soft Ultrasound. The application of soft Ultrasound is absolutely harmless for humans and animals, irrespective of age. Research has proven that ultrasound with a maximum output of 50 watt is harmless for humans and animals. The Ultrasound of the emmi®-pet uses only 0.2 Watt.

The positive effect and safety have been proven and confirmed by veterinary dentists and veterinarians.

The dental and mouth care with Ultrasound have been applied for years, without any side-effects.

Millions of Ultrasound waves are generated every minute by the patented Ultrasound chip in the brush head. The ultrasound waves are transmitted via each

bristle to the Ultrasound toothpaste. The supplied emmi®-pet ultrasonic toothpaste contains microscopic air bubbles (microbubbles) that have been enriched in an additional production step. The Goethe University of Frankfurt am Main has confirmed this. More than 3.5 million microbubbles were detected in a range of one mm³. The ultrasonic vibrations destroy the structure of the microbubbles and they collapse. This creates a suction effect that removes the biofilm and the plaque from the teeth. Afterwards we recommend brushing the bristles over the teeth to remove the loosened biofilm and the plaque from the tooth. All of this happens without your animal noticing it. The emmi®-pet is an Ultrasound unit for dental and mouth care in animals on the basis of the emmi®-dent Professional, which since 2009 has been supplied to hundreds of thousands satisfied customers worldwide.

The emmi®-pet ultrasound removes bacteria on and deep in the gums

100% ULTRASOUND TO PREVENT DENTAL DISEASES

emmi®-pet ZA Ultrasonic Toothpaste for Animals

emmi®-pet ZA - Ultrasonic toothpaste especially developed for the needs of animals.

emmi®-pet A1/A2 Ultrasonic Attachments

Large Attachment for big dogs and horses.

Small Attachment for small dogs and cats.

On hygienic grounds, a same attachment must not be used for several animals.

Replace the attachments every 3 months, because the Ultrasound softens the

bristles after this period. When the bristles are softer, the Ultrasound transmission is reduced making the tooth cleansing less efficient. Spare attachments are available in the online shop www.emmi-pet.com

Battery Care and Disposal

The battery is contained in the watertight handpiece of the unit. After use, the handpiece is to be repositioned on the charging station. When it is fully charged, it can be used for approx. 10 days without recharging. If the complete daily cleansing time exceeds 10 minutes, recharge the batteries within a shorter period. Make sure that the batteries are never completely uncharged. In case the batteries have reached a critical level, when the unit is switched on, the green display will not light up but start blinking. In such case, recharge the unit immediately in order to avoid damages to the battery. Although the unit contains a special environment-friendly lithium battery, do not dispose it with the normal household waste. For disposal, community collection points as well as other containers are available (e.g. at all sale points for similar batteries).

Warranty

For this product applies a 2 year warranty from the date of purchase. Defects which occur during the warranty period, either material and/or production defects will be repaired or the device will be replaced at no costs. As far as legally permissible, all other claims for compensation are excluded. The warranty doesn't apply in the event of intervention by third parties or disassembly by third parties without prior written approval by Emmi Ultrasonic GmbH. The warranty also doesn't apply in the event of arbitrary damage or incorrect handling. In the event of a delivered warranty adjustment, the warranty period doesn't extent. Emmi Ultrasonic GmbH reserves the right to invoice costs, incurred for unjustified complaints such as wrong installation or operating errors. For making your warranty claim, it is necessary to store and provide the receipt/proof of purchase stating the date of purchase. For a quick processing of your warranty claim please send us the device, a copy of the invoice and a short written description of the defect. Please send your claim sufficiently prepaid directly to Emmi Ultrasonic GmbH. Non-prepaid shipments can't be processed. Herewith we point out the legal rights of the consumer as well as that these are not limited by the warranty.

Thank you for your cooperation and understanding.

Emmi Ultrasonic GmbH

Gerauer Strasse 34, D-64546 Moerfelden-Walldorf

Tel.: +49 (0) 6105 - 40 67 94

info@emmi-pet.de

Selling Date

Stamp of Seller

Modifications are subject to change

© Copyright EMAG AG. All rights reserved. Made in Germany/EU. The content of this user manual is subject to the copyright. It's forbidden to copy the manual in parts and/or in total. All changes, cutbacks, enlargements or amendments, publications or translations by third parties are subject to prior written approval of EMAG AG.

Editorial status: November 2022 LK

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© Copyright EMAG AG. Alle Rechte vorbehalten. Hergestellt in Deutschland/EU. Der Inhalt dieser Bedienungsanleitung unterliegt dem Urheberrecht. Die Bedienungsanleitung darf weder ganz noch teilweise kopiert werden. Veränderungen, Kürzungen, Erweiterungen und Ergänzungen sowie jede Veröffentlichung oder Übersetzung durch Dritte bedarf der vorherigen schriftlichen Einwilligung durch die EMAG AG.