emmi[®] SKIN

SPORT

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Users' Guide

CORRECT PLACEMENT

**PATENTED AND TESTED TECHNOLOGY
MADE IN GERMANY**

emmi[®]-skin SPORT – 100% Ultrasonic for Young, Beautiful, Healthy Skin

Congratulations! You have chosen a high-quality, innovative product. In order to achieve the best possible results with your new emmi[®]-skin Ultrasonic Unit, please read the following users' guide carefully.

The emmi[®]-skin SPORT Set includes:

1 emmi[®]-skin SPORT hand unit, 1 Size L ultrasonic cover, 1 charging station, 1 emmi[®]bioMYO (optional), 1 The Thermo-effect foam (optional)

Before Turning the Unit On

Plug the charging station into a plug socket. Place the hand unit on the charging station and charge it for roughly 24 hours before using it for the first time. If it is used daily for roughly 6 minutes, the battery charge will last roughly 8 days. Afterwards, re-charge the unit for 24 hours. The green LED will light up when the charging station connection has been made properly.

Before Initial Usage

Two additional green LEDs on the hand unit indicate the battery charging state. If your emmi[®]-skin SPORT ultrasonic unit has been charged properly, place the cover perpendicularly over the hand unit from above. When removing the ultrasonic cover, pull it off perpendicularly as well.

Activation and Deactivation

Activation:

Activate the unit by pressing the power button (for roughly 1 second). Optionally, you can deactivate the tone indicating a change in position by pressing the power button for roughly 9 seconds. The green operational indicator LED will light. The battery charging state will be indicated by two other blinking green LEDs. The unit will also vibrate as an additional functional indicator for the user.

Deactivation, Automated or Manual

Deactivate the unit by pressing the power button (for roughly 1 second). The unit will automatically turn off after 6 minutes (Power Saver mode for shipping). The unit also has an electronic system for detecting the cover. If the ultrasonic cover has not been put in place, is incorrect or defective, the unit will automatically deactivate and the lower green LED will turn off.

LED Indicators

1 Green POWER LED

Press the power button for roughly 1 second (activation protection) to turn the unit on. The lowest green LED will light (operational indicator).

2 Green Fully Charged LED

This blinking green LED indicates that the battery has been fully charged.

3 Green Partially Charged LED

This blinking green LED indicates that the battery has been charged to the 30% level.

4 Blue Charging LED

If power is being properly supplied to the hand unit at 100-240V, this LED will be blue. Once the unit has been completely charged, the LED will be darker (charging state check).

5 CHECK/FAILURE Red LED

If the battery is empty or there is a failure in the unit, this LED will be red. When the unit turns off, the red LED will blink for 5 seconds and then turn off.

Using the emmi[®]-skin SPORT with 100% Original emmi[®]-skin Ultrasonic

When the unit has been turned on, you will feel it vibrate and hear a soft indicator tone. Because ultrasound cannot be seen, heard or felt, a small motor has been built into the hand unit so that the user can feel and hear when the unit has been activated and is ready for use. A soft tone will be generated every 6 seconds over the ideal usage period as well as when its position is changed. The motor will also stop briefly so that you receive tactile feedback. The tone can be turned off, as described above.

Usage

After leisure and sport activities for regeneration or in case of injuries respectively tension, we recommend you to use the emmi[®]-skin SPORT 100% Original Ultrasonic.

- 1 For reasons of hygiene, we recommend cleaning the ultrasonic cover before each use. For this purpose, we recommend using our disinfecting spray from the Emmi[®]-Beauty line, any other disinfecting spray or simple water.
- 2 Apply on the skin areas to be treated, our emmi[®]-skin crème products, as emmi[®]-bioMYO and emmi[®]-Thermo-foam. The effect of the Original 100% emmi[®]-skin Ultrasonic unit is only guaranteed in combination with the use of emmi[®]-skin crème. You will find more information about the Ultrasonic Crème from emmi[®]-skin in the emmi[®]-skin Ultrasonic Crèmes & Foams section of this users' guide or at www.emmi-skin.de.
- 3 Turn the unit on and hold it against the location for at least 6 seconds without moving it. A longer application period will be harmless, because the emmi[®]-skin SPORT uses a soft level of ultrasonic. We recommend an application period till 30 minutes- medium application duration aprox. 15 minutes- for each primary area and problem. Longer application can lead to faster, better results.
- 4 Move to the next primary skin area and repeat Step 3.
- 5 Once the primary cleaning process has been completed, the emmi skin SPORT unit will turn off automatically. Clean the ultrasonic cover and the hand unit carefully with some water or a disinfecting spray. Afterwards, you can return the emmi[®]-skin SPORT unit to the charging station with the cover on the unit.

Effect

The emmi[®]-skin SPORT unit has been developed by German engineers according to the most recent scientific discoveries. Unlike other traditional skin care devices, emmi[®]-skin SPORT does not work by means of mechanical friction or abrasively like peeling. Instead, it uses the most modern ultrasonic technology, which does not abrade the skin. This technology has been used with wonderful success in many medical areas. Using ultrasound for skin care is absolutely harmless to people. The fact that using ultrasound at up to 50 watts is harmless to the human body has been proven scientifically. The emmi[®]-skin SPORT uses 0.2 watts.

Dermatologists have tested the effects and safety of the technology. The patented piezo chip in the cover generates the ultrasonic waves and cause millions of micro-bubbles in the emmi[®]- Ultrasonic crème, which are channelled deep into layers of the skin. Human skin consists of 3 functional layers:

- the epidermis
- the dermis
- the hypodermis

Conventional crèmes only affect the epidermis. In contrast, the emmi[®]-skin SPORT Ultrasonic unit successfully channels the contents of emmi[®]-skin crèmes into the hypodermis and the leather of the skin due to its unique effect. Doing this provides the skin with effective support as it forms new cells and regenerates itself. The positive effect of using ultrasound also releases the body's own energies, improves the transportation mechanisms of the lymph nodes and improves circulation.

People are often told that water is required as a transfer medium. The reality is that the appropriate transfer medium depends on the frequency and performance rate of the ultrasonic waves. Medical applications are well aware of these issues and therefore use salves and gels with different consistencies as appropriate transfer media. We have developed a special ultrasonic crème for injuries and tense muscles. To achieve the ideal results, please use emmi[®]- skin Crème together with emmi[®]-skin SPORT.

emmi[®]-skin Ultrasonic Cremes & Foams

Thermo-effect foam:

The Thermo-effect foam acts within seconds as heat therapy for joint and muscle pain. The warmth simulates the blood circulation, which proves itself as very successful to muscular pain. Through the heating effect the pain can be relieved and the muscles are warmed and relaxed (please follow the instructions use of the product).

emmi[®] bioMYO:

The biologically active components of hellebore, horse chestnut and juniper extracts have a relaxing, anti-inflammatory, decongestant effect and promote capillary blood circulation. Due to the stimulating effect of blood circulation, emmi[®] bioMYO is suitable for a series of pain disorders like arthroses, myalgia, pulled muscle, swelling and other various neuralgias. Please note that emmi[®]-skin optimised products will also work without the emmi[®]-SPORT ultrasonic unit. Of course, the emmi[®]-SPORT unit will increase the positive effect of the crèmes.

Replacement Covers

You should be aware that the use of the ultrasonic cover by multiple people may not be appropriate for reasons of hygiene and it should be replaced after roughly three months of use. Only one person should use a given ultrasonic cover. You can purchase replacement covers at the same place where you purchased your emmi[®]-skin SPORT unit.

Battery Charging Technology

The hand unit can be put back on the charging station after use to re-charge the lithium battery or to be able to use the-emmi[®]-skin SPORT unit with a full battery charge every 8 days. If the period of usage is more than 6 minutes per day, the battery should be recharged for a brief period. If the battery is almost empty, the operational indicator (lowest green LED) will not light up. Instead, the red LED will light up, blink for 5 seconds and the unit will be completely turned off. If that happens, the unit must be recharged for 24 hours.

Replacing the Battery & Protecting the Environment

The unit has a special, long-life, environmentally friendly lithium battery located in the hand unit. If the red LED continues to light up after long periods of use, despite recharging, the battery needs to be replaced. The battery can only be replaced by an emmi[®]- skin SPORT lithium power-saver battery (original replacement part with an internal fuse). We are able to replace the battery for you, if you do not want to

do this by yourself. If you want to replace the battery yourself, follow these instructions.

- 1 Loosen the screw at the bottom of the hand unit and pull the battery holder out with the floor plate.
- 2 When inserting the new battery, check the polarity. Match the battery's positive terminal with the floor plate holder.
- 3 Be sure that all of the seals sit properly when re-assembling the parts. Doing this ensure the IP54 class protection. You can purchase replacement batteries at the same place where you purchased your emmi®-skin SPORT unit.

Note: any guarantee claims will only be applicable if the product has been treated properly and original parts have been used when replacing the battery.

Disposal of Batteries

Do not put the special battery in the unit in the household garbage. To dispose of the battery, please use the mutual return collection containers or other return centre that are available (such as any point of purchase for similar batteries). Please note the dangers indicated below.

Dangers

The emmi®-skin SPORT unit includes a battery with a built-in fuse. Despite this, do not give the battery to children and comply with the notes about disposal above. If liquids leak from the battery despite the safety measures taken, do not allow them to contact the skin and dispose of the battery as described above. The connection with wall socket also requires safety measures. The charging station has been cast to insulate against electricity, however the charging station should never be placed in water. Make sure that the power cable for the charging station is free of damage. If the housing, plug or power cable becomes damaged, the entire charging station should not be used anymore and needs to be replaced. The emmi®-skin SPORT unit is only appropriate for the usage described above.

Guarantee

The guarantee will only apply to the 24 months after the date when the emmi®-skin SPORT was purchased; the battery guarantee only applies for 6 months after the date of purchase. Defects in the form of material or manufacturing defects that appear during the guarantee period will be corrected free of charge. The correction may take the form of repairs or replacement. To the extent legally allowed, all other claims for compensation shall be excluded. The guarantee will be voided if the unit is opened by a third party or disassembled without the prior permission of Emmi Ultrasonic GmbH.

The guarantee will also be voided by intentional damage or incorrect usage. The guarantee period does not imply the provision of a guarantee. We retain the right to bill for any expenses that arise from unjustified complaints, such as errors made during installation or operation. The purchase receipt with the date of purchase must be kept to make guarantee claims. So that we can process your guarantee claim quickly, please the entire unit, including the ultrasonic cover and the charging station with sufficient postage directly to Emmi Ultrasonic GmbH. Shipments that are not free of charges will not be accepted. Thank you for your understanding.

Proper Usage

The products may only be used for the issues intended for in the catalogues and according to the technical descriptions, meaning cleaning and skin care purposes, and only in connection with the equipment and components recommended and allowed by Emmi Ultrasonic GmbH. ⚠ Any usage other than this or beyond this shall be considered as unintended. Emmi Ultrasonic GmbH will not assume any liability for improper usage, any damage arising and does not guarantee the problem-free function of the components. Compliance with the users' guide is a component of proper usage. The symbol on the name plate also refers to this. The smooth and safe operation of the products assumes proper transportation, storage, placement and assembly, as well as careful operation and maintenance. Compliance with repetitious checks in accordance with BGV A2 is mandatory. The unit must be setup so that the charging station's plug will be freely accessible at all times.

Improper Usage

Danger to life and limb as well as material damages may be caused by improper usage of the product. Improper usage is forbidden for that reason. Emmi Ultrasonic GmbH will not assume any liability for improper usage, any damage arising and does not guarantee the problem-free function of the components. Compliance with the following instructions is mandatory to exclude improper usage of the unit:

- The emmi[®]-skin SPORT unit for skin cleaning using Original 100% emmi[®]- skin Ultrasonic may only be operated using the voltage that has been printed on the name plate (on the bottom of the charging station).
- Do not use combustible or aggressive chemicals, acids and similar agents for cleaning.
- If damage occurs during transportation do not connect the charging station to the plug socket.
- Only use the unit in internal rooms at room temperature.

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Emmi Ultrasonic GmbH • Gerauer Straße 34